

2024-04-26, v2.1.0

1. #34635 Add BeaconNext to EVAPlayer

2023-03-15, v2.0.8

1. #25171 New parameter for ENTER key in scanOnlyMode

2022-09-23, v2.0.7

1. #21232 MediaProjectionManager change on A10 and A11

2022-08-17, v2.0.6

1. #20340 Not pairing Pro Glove Intent Scanner

2022-05-23, v2.0.5

1. #14280 Android 10 - 11
2. #18083 ProGlove scanner configuration problem

2021-05-19, v2.0.4

1. #1747 Better startup management with EVA Selector
2. #1761 Add new ViaIntent solution for Unitech scanner

2021-04-08, v2.0.1 - v2.0.3

1. #1533 Prise en compte de la Bague ProGlove - it didn't work with some configurations.
2. #1733 Add support for IntentScanner
3. #1730 Android 10 -11 (partly fixed, changed permission for reading /sdcard)

2021-03-16, v2.0.0

1. #1717 New APK certificate
2. #1533 Prise en compte de la Bague ProGlove

2020-05-27, v1.4.1

1. #771 Adapt Eva Selector/Player/Logger to new API 29 (Android 10) about app launch

2019-10-29, v1.4.0

1. #854 - Create a text file with detailed version when launching

2019-05-31, v1.3.10

1. Fixed bug, when voixtreme folder was missing on device. Did not notify about file copy error.
2. Bug fix, If you install a Bar scanner app without starting it first once, sometimes it does not respond to a BroadCast. Warning dialog added for the user.

2019-05-16, v1.3.9

1. Version name change for installation script.

2019-05-06, v1.3.8.3

1. When "run" application was missing, Selector did not close down automatically

2019-04-30, v1.3.8.2

1. Fix bug when opening EVA Selector without having the proper configuration on the device.
2. Fix bug when opening Calibrate, the vocalized APP opened just after.

2019-04-29, v1.3.8.1

1. Incremented the delay when opening the APP to vocalize